

NW 200 nutch filter unit



Characteristics:

- batch-wise filtration of pharmaceutical solutions
- parts in contact with product made entirely of 3.3 borosilicate glass and PTFE
- 15 l batch volume
- DN 200 filter plate and vessel
- optional suction operation at excess pressure, design with splinter shield and perforated guard plate with good visibility
- manual agitation with process-specific height-adjustable stirrer blade
- dairy couplings for feed and solvent nozzles, with blank cover
- pressure, pH and temperature measurement
- bottom outlet valve
- combined feeder tube and baffle
- mobile stainless steel frame for easy assembly and disassembly, manual lowering of filter vessel, tilting construction,
- completely preassembled package unit
- directives
 - pressure equipment directive 97/23 EC
 - FDA material certificates