

## 200 l autoclavable mixing vessel, pharma design



mixing vessel 200 l



complete unit in autoclave

### Characteristics:

- mixing vessel with
  - 200 l spherical vessel, 3.3 borosilicate glass
  - agitator with single-acting mechanical seal, manually adjustable rotation speed 0-700 rpm
  - internal agitator bearings, corrosion-proof at product contact
  - 3-stage propeller stirrer, PTFE at product contact
  - 3.3 borosilicate glass inlet tube for gas dispersal
  - **NORMAG** angle valves in autoclavable, valve upper part in stainless steel design
  - low-dead-space bottom outlet valve in autoclavable, valve upper part in stainless steel design
  - TriClamp PTFE-Adapter for quick and easy connection by the client
  - DN 150 quick release
- mobile stainless-steel frame for rail systems, suitable for the client's own autoclaves
- agitator motor easy to install and remove with wall-mounted lifting tackle
- unit without motor completely autoclavable
- all flange connections in stainless steel
- directives
  - CE compliance / machinery directive
  - pressure equipment directive 97/23 EC
  - FDA material certificates
  - class B clean rooms