

200 l stirring vessel with insertion heat exchanger



Characteristics:

- 200 l spherical vessel with stirrer for mixing and tempering
 - no heating jacket
 - floor drain valve with limited dead space, fully drainable
 - for pharmaceutical and chemical uses
- Stainless steel frame
- Tempering via internal spiral heat exchanger with
 - internal forced re-circulation and correspondingly good exchange of substances
 - increased heat transfer
 - with higher permitted temperature range
- stirrer with PTFE jacketed impeller stirrer for internal forced re-circulation of liquid in the 200 l sphere, single-acting, dry running, FDA-compliant and corrosion-resistant GLRD, up to 700 rpm, manual speed setting
- Stainless steel frame, mobile
- Optional tempered discharge
- Optional insertion heat exchanger in stainless steel, Hastelloy or Tantal
- Directives
 - Machinery Directive/CE conformity
 - Pressure Equipment Directive 97/23 EC
 - FDA Material Certificate
 - TA Air Certificate